

PerSeptive Biosystems  
Voyager Elite

MALDI-TOF  
Mass  
Spectrometer





Figure 1-5 Voyager-Elite Mass Spectrometer

**General Information file**  
*(do NOT remove)*

CLEAR TOPPER

Parts of the system include:

- **Mass spectrometer**—A time-of-flight mass spectrometer, described in Section 1.3.2, Mass Spectrometer. The high-vacuum system of the mass spectrometer is described in Section 1.3.3, Vacuum System.
- **Computer/Data System**—A computer that operates the Voyager control software and the Voyager processing software. You control the mass spectrometer using the computer.
- **Digitizer**—An analog-to-digital converter that allows the signal from the mass spectrometer to be transferred to the computer.
- **Video Monitor**—A monitor that displays a real-time video image of the sample spot for examination during acquisition.
- **Control Stick**—A peripheral device that controls position of the sample plate in the mass spectrometer. Allows you to start and stop acquisition and transfer data to processing software.



- = V to MS
- - - = laser beam
- · - = addr/data bus
- ▲— = trigger out

CABLE DISTRIB.

CONTROL BOARD

TO PC

to Digitizer

LASER

HF Counter

HF Counter

Detector Delay Timer

30 kv SRC HV

HV SW (DE)

30 kv Grid HV

+/- 100 V G. Wire

HV SW & 2 kv

Mirr. HV Driver

Mirr. HV

30 kv Mirr. HV

Lin. HV Driver


Lin. HV

DAC

DAC

DAC

DAC



**Figure 1-3 Time of Flight Analysis**

The time required for ions to reach the detector at the opposite end of the flight tube is measured. The number of ions reaching the detector at any given time is also measured, and is referred to as ion intensity (abundance) or signal intensity.

Drift time is proportional to the square root of the mass and is calculated by the following equation:

$$t = s \left( \frac{m}{2KE} \right)^{1/2}$$

where:

- t = drift time
- s = drift distance
- m = mass
- KE = kinetic energy

$$E_c = \frac{1}{2} m v^2 = \frac{1}{2} m \left( \frac{s}{t} \right)^2$$

$$\Rightarrow t = s \left( \frac{m}{2E_c} \right)^{1/2}$$

# MATRIX-ASSISTED LASER DESORPTION IONIZATION-TIME OF FLIGHT (MALDI-TOF) ANALYSIS OF PEPTIDES USING POST-SOURCE DECAY

A.M. Fallick, PerSeptive Biosystems Western Technical Center, San Francisco, CA


*In recent years, mass spectrometry has become an essential tool for the analysis of biopolymers (1-3). Accurate molecular weight information can be determined in a few minutes on extremely small quantities of material in a straightforward fashion. In addition, it is also possible to obtain sequence information using various mass spectrometric methods such as MS/MS (mass spectrometry/mass spectrometry). This Technical Bulletin describes the use of post-source decay for peptide analysis.*

## MS/MS ANALYSIS

MS/MS analysis requires the use of two stages of mass separation. A four-sector mass spectrometer consisting of two double-focusing magnetic sector instruments placed end-to-end is frequently used for this type of analysis (4,5). The first instrument is used to select the ions to be analyzed; in the case of a peptide mixture,

ions formed from a selected peptide are allowed to pass into the second mass spectrometer, while ions of the "wrong" mass are rejected. A collision cell filled with a low pressure of an inert gas such as helium is placed between the two mass spectrometers. Selected ions pass into the cell and collide with gas molecules; as a result some ions decompose due to the internal energy gained in the collision. The fragments that are formed are mass-analyzed in the second mass spectrometer. The resulting mass spectrum, often called a collision-induced dissociation (CID) spectrum, yields a great deal of information about the starting compound.

MS/MS analysis can also be performed using quadrupole instruments referred to as triple quadrupoles (6). The first quadrupole selects the ion of interest, the second is a collision cell, and the third is used to analyze the fragments. Quadrupoles generally operate



**Figure 1.** Principle of a linear TOF mass spectrometer. Ions are formed at the target by a laser pulse striking a mixture of sample and matrix placed on the sample plate (target).

at lower energies than sector instruments and while the fragmentation is often simpler, it may be less informative.

## POST-SOURCE DECAY (PSD) ANALYSIS


The post-source decay method (7) used in time-of-flight (TOF) mass spectrometers has much in common with MS/MS in that fragments of a selected analyte molecule are formed and analyzed in the instrument. This method has been successfully used to determine and verify peptide sequences (8,9); the spectra obtained are very similar to those obtained using the double mass spectrometer techniques (10).

In traditional mass spectrometric analysis, analyte ions are normally produced by matrix-assisted laser desorption ionization (MALDI) (11, 12) yielding very high sensitivity results. The basic TOF instrument (without PSD) operates as shown in Figure 1. Ions are formed in the ion source by a laser pulse, then rapidly accelerated to a high kinetic energy (20-30 keV). In a linear TOF instrument, once ions leave the short acceleration region, they travel down a 1-2 meter tube to the detector, where their arrival is timed. Mass-to-charge ratio,  $m/z$ , is determined from the time elapsed from ion formation (laser pulse) to ion arrival time at the detector. Ions typically acquire enough internal energy at the time of ionization to undergo unimolecular


decomposition as they travel along the flight tube; in a linear instrument, however, both intact ions and any fragment ions and neutrals arrive at the detector at virtually the same time. As a result, only one peak is detected, at a time corresponding to the precursor ion mass.

Some TOF instruments are equipped with an electrostatic mirror (reflector) (13) that reflects ions back towards a second detector (Figure 2). The original function of the reflector is to reduce the translational kinetic energy spread of the ions so as to improve resolution. Faster ions penetrate farther into the reflecting field and thus travel a longer distance than slower ions. As a result, if instrument parameters are set properly, all ions of the same  $m/z$  arrive at the detector at nearly the same time.

It is important to realize that all ions formed in the source and accelerated normally will have approximately the same amount of translational kinetic energy



**Figure 2.** Diagram of a reflecting type TOF mass spectrometer. This instrument is equipped with a beam guide for more efficient ion transmission and a timed ion selector.



**Figure 3.** Electrostatic mirror (reflector) portion of a TOF instrument. An ion fragments before entering the mirror (anywhere in the field free region) and the charged and neutral fragments enter the mirror together. The neutral fragment strikes the detector at the back of the mirror (normally turned off in this mode) and is lost. The charged fragment turns around in the mirror and is detected at the detector located at the front of the mirror.

and thus penetrate about the same distance into the mirror. Ion speeds however, differ according to their  $m/z$  values and therefore, so do their arrival times. This includes fragment ions formed concurrent with or immediately after ionization ("prompt" fragments).

Fragment ions formed as a result of unimolecular decay after acceleration of the precursor ion out of the ion source and before reflection (post-source decay), are quite different. These ions do not follow the same path as fully accelerated ions and as a result, reach the reflector detector at a different time than either their intact precursor ions or normal fragment ions formed in the source. Neutral fragments are not reflected and are not detected at the reflector detector (Figure 3).

Consider the example of a singly charged ion of mass 1000 that decomposes to an  $m/z$  500 ion in flight after it leaves the ion source. The precursor ion ( $m/z$  1000) will reach the reflector detector at the expected time for its mass; the fragment ion will have less kinetic energy (half as much in this example) and will be moving more slowly and thus not penetrate as far into the mirror as would a fully accelerated ion of the same mass. While these properties (slower speed, shorter path) have opposite effects on the arrival time, the result is that the ion arrives at the detector at a slightly shorter time than a normal fragment of the same mass. A fragment

of  $m/z$  250 from the same precursor would have even less kinetic energy and would penetrate the mirror even less than the  $m/z$  500 ion. However, given the instrument parameters, one can accurately calculate the mass of a fragment ion formed in flight from its arrival time. The resulting spectrum contains peaks at the masses of fragments formed by PSD.

Ions that do not penetrate fully into the mirror are not refocused as effectively as ions that travel almost the entire mirror length. Ions that have less kinetic energy than this optimum value can still be correctly focused simply by making the reflecting field weaker. The kinetic energies of PSD fragment ions vary according to their masses, and as a result, no one reflector field will be optimal for all ions. The range of energies that can be properly focused is such that if a mirror voltage is varied in 5-10 steps, all of the ions of interest can be satisfactorily focused and detected. Reducing the mirror voltage has an additional advantage: Fully accelerated ions are not reflected (they pass through the back of the mirror) and as a result, do not appear in the spectrum; this includes prompt fragment ions formed in the source and non-fragmenting ions from other species that may be present. The PSD spectrum is formed by assembly of the 5-10 steps, normally accomplished by the data system of the mass spectrometer.




Figure 4. PSD spectrum of ~2 picomoles of the hemoglobin tryptic peptide VHLTPEEK. Peaks in the spectrum are marked in accordance with the standard nomenclature for fragment ions (14).

A TOF instrument that is properly equipped to do PSD analysis usually includes a timed ion selector (Figure 2). This device, located along the flight path of the precursor ions, allows the selection of precursor ions of a certain mass so that the PSD spectrum contains only fragments from the selected precursor ion. As a result, a pure PSD spectrum of one component of a mixture can be obtained without having to perform a prior separation. The timed ion selector works simply by opening an electrostatic gate at the arrival time of the desired precursor ion. These ions pass through normally; ions arriving at different times (different  $m/z$  values), however, are not allowed to pass. Typically, the gate can be set to discriminate against ions that differ from the desired ion by approximately one percent or more.

Peptides typically fragment to give both sequence ions and amino acid specific ions at low masses. Figure 4 shows the PSD spectrum of a hemoglobin tryptic peptide (beta T1). The spectrum shows a complete series of  $y$  ions (cleavage between the carbonyl and amino groups in the peptide backbone with charge retention on the C-terminal fragment). From differences in mass between the members of this series of ions, one can determine or verify the amino acid sequence of the peptide. The low mass ions marked with single letter codes are characteristic of the presence of specific amino acid residues in the peptide and provide information as to the amino acid composition of the peptide.

## REFERENCES

1. Biemann, K. 1992. Annual Review of Biochemistry **61**:977-1010.
2. Burlingame, A.L., R.K. Boyd, S.J. Gaskell. 1994. Analytical Chemistry **66**:634R-683R.
3. Chiat, B.T., S. Kent. 1992. Science **257**:1885-1894.
4. Walls, F.C., M.A. Baldwin, A.M. Falick, B.W. Gibson, S. Kaur, D.A. Maltby, B.L. Gillece-Castro, K.F. Medzihradsky, S. Evans, A.L. Burlingame. 1990. in Biological Mass Spectrometry. Ed., Burlingame, A.L., J.A. McClosky (Elsevier, Amsterdam) pp. 129-146
5. Ashcroft, A.E., P.J. Derrick. 1991. in Four-sector tandem mass spectrometry of peptides. Ed., Desiderio, D.M. (CRC Press, Boca Raton FL) pp. 121-138.
6. Yost, R.A., R.K. Boyd. 1990. Meth. Enzymology **193**:154-200.
7. Spengler, B., D. Kirsch, R. Kaufmann, E. Jaeger. 1992. Rapid Communications in Mass Spectrometry **6**:105-108.
8. Kaufmann, R., B. Spengler, R. Lutzenkirchen. 1993. Rapid Communications in Mass Spectrometry **7**:902-910.
9. Kaufmann, R., D. Kirsch, B. Spengler. 1994. Int. J. Mass Spectrom. Ion Processes **131**:355-385.
10. Rouse, J.C., W. Yu, S. Martin. 1995. J. Am. Soc. Mass Spectrometry **6**:822-835.
11. Karas, M., D. Bachmann, U. Bahr, F. Hillenkamp. 1987. Int. J. Mass Spectrom. Ion Processes **78**:53-68.
12. Hillenkamp, F., M. Karas, R. Beavis, B.T. Chait. 1991. Analytical Chemistry **63**:1193A-1203A.
13. Mamyrin, B.A., V.I. Karataev, D.V. Schmikk, V.A. Zagulin. 1973. Sov. Phys. JETP **37**:45.
14. Biemann, K. 1990. Meth. Enzymology **193**:886-887.



## PerSeptive Biosystems

### U.S. HEADQUARTERS:

PerSeptive Biosystems, Inc.  
500 Old Connecticut Path  
Framingham, MA 01701  
Telephone: (508) 383-7700  
Toll-free: (800) 899-5858  
Fax: (508) 383-7881  
<http://www.pbio.com/>

### JAPAN:

Nihon PerSeptive Ltd.  
5 Koike Building  
1-3-12, Kita-shinagawa 1-chome  
Shinagawa-ku,  
Tokyo, 140 Japan  
Telephone: +81 33 471 8191  
Fax: +81 33 471 7035

### EUROPEAN HEADQUARTERS

(GERMANY):  
PerSeptive Biosystems GmbH  
Otto-van-Guericke-Ring 7  
D-65205 Weisbaden  
Germany  
Telephone: +49 6122 9830  
Fax: +49 6122 9830 30

### UNITED KINGDOM:

PerSeptive Biosystems UK Ltd.  
3 Harford Court  
Foxholes Business Park  
John Tate Road  
Hertford SG13 7NW  
Telephone: +44 1992 507100  
Fax: +44 1992 553858

### FRANCE:

PerSeptive Biosystems France Ltd.  
164 avenue Joseph Kessel  
Parc d'activité "Parkile II"  
F-78960 Voisins le Bretonneux  
Telephone: +33 1 34523030  
Fax: +33 1 30640909

Perfusion Chromatography and POROS are registered trademarks of PerSeptive Biosystems, Inc. BioCAD, ID, INTEGRAL, ImmunoDetection, ConSep LC100, MemSep, Poroszyme, Real-Time Process Monitor, Voyager, Biospectrometry, and 9050 Plus PepSynthesizer are trademarks of PerSeptive Biosystems, Inc.





POROS® products and perfusive chromatography are covered by US patents 5,000,352; 5,019,270; 5,228,989 issued to PerSeptive Biosystems, Inc.; foreign patents pending. Chromatography materials of certain pore geometries, including wide pore supports, can permit perfusive chromatography, that is, separations where the rate of intraparticle convection exceeds the rate of diffusion. Use of any such chromatography media that permit performance at high linear velocity without license, constitutes infringement of US patent 5,019,270 and US patent 5,228,989. A limited license to use the patented perfusive chromatography process is granted with the purchase of POROS® products from PerSeptive Biosystems. The license may not be transferred except with the transfer of POROS® columns or media, and terminates upon expiration of the useful life of the product.






http://matrix.depot

FOR ADDITIONAL INFO. TRY THE WEB:

| Matrix   | Application  | Structure       | Color of crystals/solution | Matrix Solution  |
|--|--|-----------------|----------------------------|--|
| Sinapinic acid<br>(see mass spectrum on page C-2)<br>MW 224.07                                       | <ul style="list-style-type: none"> <li>• Peptides</li> <li>• Proteins</li> </ul>   | <p>PBI00251</p> | White                      | <ul style="list-style-type: none"> <li>• 10 g/liter in 50:50 water/ethanol<sup>1</sup></li> <li>• 10 g/liter in 70:30 water/acetonitrile (0.1% TFA)<sup>2</sup></li> </ul>     |
| Alpha-cyano-4-hydroxy cinnamic acid ( $\alpha$ CHCA)<br>(see mass spectrum on page C-2)<br>MW 189.04 | <ul style="list-style-type: none"> <li>• Peptides</li> <li>• Proteins</li> </ul>   | <p>PBI00252</p> | Yellow                     | <ul style="list-style-type: none"> <li>• 10 g/liter in 70:30 water/acetonitrile (0.1% TFA)<sup>3</sup></li> <li>• 10 g/liter in 50:50 water/acetonitrile (0.1% TFA)</li> </ul> |
| 2,5-dihydroxybenzoic acid (2,5-DHB)<br>MW 154.03   | <ul style="list-style-type: none"> <li>• Peptides (mixtures)</li> <li>• Carbohydrates</li> <li>• Glycolipids (neg ion mode)</li> <li>• Polar synthetic polymers</li> </ul> | <p>PBI00253</p> | White                      | 10 g/liter in 90:10 water/ethanol <sup>4</sup>   |

| Matrix   | Application   | Structure   | Color of crystals/solution | Matrix Solution  |
|--|---|---|----------------------------|--|
| 2-(4-hydroxy-phenylazo)-benzoic acid (HABA)<br>MW 242.07 | <ul style="list-style-type: none"> <li>• Proteins</li> <li>• Polar and nonpolar synthetic polymers</li> </ul> |  <p style="text-align: right;">PB100255</p>   | Orange                     | ~1.3 g/liter in 50:50 water/acetonitrile or in 40:40:20 water/acetonitrile/methanol <sup>5</sup> |
| 2-aminobenzoic acid (anthranilic) acid<br>MW 137.05      | Oligonucleotides  |  <p style="text-align: right;">PB100256</p>   | Yellow                     | ~10 g/liter in 80:20 water/acetonitrile (20% w/w nicotinic acid can be added) <sup>6</sup>       |
| 3-hydroxypicolinic acid<br>MW 139.03                     | Oligonucleotides  |  <p style="text-align: right;">PB100257</p>  | Light brown                | ~80 g/liter in 50:50 water/acetonitrile <sup>7</sup>   |
| Dithranol<br>MW 226.06                                   | Nonpolar synthetic polymers<br>+ polar polymers   |  <p style="text-align: right;">PB100258</p> | Yellow                     | 10 g/liter in tetrahydrofuran <sup>8</sup>   |

| Matrix   | Application      | Structure   | Color of crystals/<br>solution | Matrix Solution  |
|--|------------------|---|--------------------------------|--|
| 2,4,6 Trihydroxy acetophenone (THAP)<br><br>MW 168.04    | Oligonucleotides |   | White                          | 0.5 M 2,4,6-THAP in ethanol, mix with 0.1 M aqueous diammonium hydrogen citrate (2:1) <sup>9</sup> |
| Picolinic acid (YAG laser only, 266 nm)<br><br>MW 123.03 | Oligonucleotides |   | White                          | 1.2 to 2.4 M (150 g/liter to 300 g/liter) in water <sup>10</sup>                                   |
| Nicotinic acid (YAG laser only, 266 nm)<br><br>MW 123.03 | Proteins         |  | White                          | 0.05M in water <sup>11</sup>   |

1. Beavis, R.C., B.T. Chait, *Rapid Commun. Mass Spectrom.* 1989, 3, 432-435
2. Beavis, R.C., B.T. Chait, *Anal.Chem.* 1990, 60, 1835-1840
3. Beavis, R.C., B.T. Chait, *Org. Mass Spectrom.* 1992, 27, 159-156
4. Stupet, K., M. Karas, F. Hillenkamp, *Int. J. Mass Spectrom. Ion Proc.* 1991, 111, 89-102
5. Juhasz, P., C.E. Costello, K. Blummann, *J. Am. Soc. Mass Spectrom.* 1993, 4, 399-409
6. Nordhoff, E. *et al.*, *Rapid Commun. Mass Spectrom.* 1992, 6, 771-776
7. Wu, K.J., A. Steding, C.H. Becker, *Rapid Commun. Mass Spectrom.* 1993, 7, 142-146
8. Juhasz, P., C.E. Costello, *Rapid Commun. Mass Spectrom.* 1989, 7, 343-351
9. Scher, M., *et al.*, *Nucleic Acids Research* 1993, 21, 3191-3198
10. Chen, C.H., *et al.*, *Rapid Commun. Mass Spectrom.* 1994, 8, 673-677
11. Karas, M., U. Bahr, F. Hillenkamp *Int. J. Mass Spectrom. Ion Proc.* 1989, 92, 231-242

Tue Jul 29 10:43:26 1997

FYI

Date: Tue, 29 Jul 1997 10:22:18 -0400  
From: HEMPEL@CPWSCA.PSC.EDU  
Subject: FYI  
To: mbier@andrew.cmu.edu  
Message-ID: <970729102218.2260928c@CPWSCA.PSC.EDU>

Hi Mark,

This came in that I thought would be of interest - you're not on the ABRF network anymore are you?

Otherwise, I do want to come over for some ESMS, but I need to get a peptide digest separated first. More on that soon.....JH

~At 9:39 AM -0800 7/28/97, Richard S. Johnson wrote:

>Dear ABRF:

> Does anyone have a favored sample preparation method for obtaining high  
>sensitivity peptide measurements on a maldi-tof?

Richard,

I've had good luck using a modified version of Xiang and Beavis' "polycrystalline" method. It seem to work best with the stipled or roughened sample wells. I simply make a 50 mg/ml solution of alphacyano in 99% or 100% acetone and spot 0.1-0.3 ul to the position. After drying (almost instantaneous) I use a piece of hard paper to crush the crystals (or just wipe with a kimwipe, personal commun. with Ron Beavis,) and then wipe as hard as I can with a folded dry kimwiper. Then I make a 1:10 dilution of the aqueous sample in saturated alphacyano in 1:2 0.1%TFA/ACN, and spot 1.0 ul. As soon as the layer appears, remove the liquid, dry, and wash with 2 X 5ul deionized water. For the highest sensitivity, I've found you can get away spotting 1.0 ul of the saturated matrix and adding no more than 0.1-0.2 ul of the sample directly to the slide. If all goes well, I generally see high resolution but like you, greatly reduced numbers of good shots per sample, i.e. 1-4 shots before the signal degrades. So, this method might not be so hot for PSD or PSD/CAD. However, with this technique and acetylated angiotensinogen, I have been able to obtain better than 10:1 S/N with <800 attomoles spotting 50-100 nanoliters.  
best of luck,

---

Gary M. Hathaway, Director  
PPMAL - Protein/Peptide Micro Analytical Laboratory  
California Institute of Technology  
139-74, Division of Biology  
Pasadena, CA 91125

<http://www.cco.caltech.edu/~ppmal>

Tue Jul 29 10:43:26 1997

FYI

email Gary: [hathaway@caltech.edu](mailto:hathaway@caltech.edu)

email facility: [ppmal@caltech.edu](mailto:ppmal@caltech.edu)

phone: lab (818) 395-6388 or office (818) 395-2769

FAX (818) 449-3414

---